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

Class	Mutation	Location of mutation	Kinase activity	Comments
1	V600E/K/D/R	Activation segment	High	Occurs most frequently in melanomas
2	R462I	G-rich loop	Intermediate	
2	I463S	G-rich loop	Intermediate	
2	G464E/V/R	G-rich loop	Intermediate	Breast cancer G464V
2	G469 A/V/S	G-rich loop	High/Intermediate	NSCLC G649 A; G649 V melanoma
2	E586K	$\beta 7$ - $\beta 8$ loop	High	
2	F595L	DFG	Intermediate	
2	L597Q/R/S/V	Activation segment	High/Intermediate	Melanoma L597R; NSCLC L597V
2	A598V	Activation segment	Active	
2	T599I	Activation segment	Intermediate	
2	K601E/N/T	Activation segment	High	Melanoma K601E
2	A727V	Distal to $\alpha 1$	Intermediate	
2	B-Raf fusions		High	
3	G466 A/E/V/R	G-rich loop	Impaired	Stomach cancer G466V
3	S467 A/E/L	G-rich loop	Impaired	
3	G469E	G-rich loop	Impaired	
3	K438M	Proximal to $\beta 1$	Impaired	
3	D594 A/E/G/H/N/V	DFG	None	Melanoma D594N
3	G596 A/C/D/R	DFG	Impaired	Colorectal cancer G596R

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






- 1 | Conde E, Hernandez S, Martinez R, et al. Assessment of a New ROS1 Immunohistochemistry Clone (SP384) for the Identification of ROS1 Rearrangements in Patients with Non-Small Cell Lung Carcinoma: the ROSING Study. *Journal of Thoracic Oncology* 2019;14:2120-2132.
- 2 | Hofman V, Rouquette I, Long-Mira E, et al. Multicenter Evaluation of a Novel ROS1 Immunohistochemistry Assay (SP384) for Detection of ROS1 Rearrangements in a Large Cohort of Lung Adenocarcinoma Patients. *J Thorac Oncol* 2019;14:1204-1212.
- 3 | Huang RSP, Smith D, Le CH, et al. Correlation of ROS1 Immunohistochemistry With ROS1 Fusion Status Determined by Fluorescence In Situ Hybridization. *Archives of Pathology & Laboratory Medicine* 2020;144:735-741.
- 4 | Davies K, T. Le A, Jamie Sheren, et al. Comparison of Molecular Testing Modalities for Detection of ROS1 Rearrangements in a Cohort of Positive Patient Samples. 2018;13:1474-1482.
- 5 | Ryma Benayed, Michael Offin, Kerry Mullaney, et al. High Yield of RNA Sequencing for Targetable Kinase Fusions in Lung Adenocarcinomas with No Mitogenic Driver Alteration Detected by DNA Sequencing and Low Tumor Mutation Burden. 2019;25:4712-4722.
- 6 | 日本肺癌学会バイオマーカー委員会. 肺癌患者における次世代シーケンサーを用いた遺伝子パネル検査の手引き, 2019年12月23日 第 1.1 版. Available at <https://www.haigan.gr.jp/uploads/files/NGS%E3%82%AC%E3%82%A4%E3%83%89%E3%83%A9%E3%82%A4%E3%83%B3%20%E7%90%86%E4%BA%8B%E4%BC%9A%E6%89%BF%E8%AA%8D%20-%201.1%281%29.pdf>.
- 7 | Roskoski R Jr. Targeting oncogenic Raf protein-serine/threonine kinases in humancancers. *Pharmacol Res.* 2018 Sep;135:239-258.
- 8 | Toon CW, Walsh MD, Chou A, et al. BRAFV600E Immunohistochemistry Facilitates Universal Screening of Colorectal Cancers for Lynch Syndrome. *The American journal of surgical pathology* 2013;37:1592-1602.
- 9 | Long GV, Wilmott JS, Capper D, et al. Immunohistochemistry is highly sensitive and specific for the detection of V600E BRAF mutation in melanoma. *The American journal of surgical pathology* 2013;37:61-65.
- 10 | Lasota J, Kowalik A, Wasag B, et al. Detection of the BRAF V600E Mutation in Colon Carcinoma. *The American journal of surgical pathology* 2014;38:1235-1241.
- 11 | Ilie M, Long E, Hofman V, et al. Diagnostic value of immunohistochemistry for the detection of the BRAFV600E mutation in primary lung adenocarcinoma Caucasian patients. *Ann Oncol* 2013;24:742-748.
- 12 | Seto K, Haneda M, Masago K, et al. Negative reactions of BRAF mutation-specific immunohistochemistry to non-V600E mutations of BRAF. *Pathology International* 2020;70:253-261.